

Interchangeable System - 100% galvanic separated (power supply and USB signal)

Video Rhino Laryngoscope



RSX-i / RSX-iP (Paediatrics)

For interchangeable Systems

Interchangeable probe RSX-i = 3.9 mm diameter

Interchangeable probe RSX-iP (Paediatrics) = 2.9 mm diameter

Connectable with: HDMI digital video processor
USB digital video processor



RSX-i	=	3.9mm CMOS 400x400 FOV 90°
RSX-iP	=	2.9mm CMOS 400x400 FOV 90°

Digital Video Processor


HDVP / UDVP

HDVP = HDMI Digital Video Processor


UDVP = USB Digital Video Processor

Input: RSX-i / RSX-iP

Output: HDMI with 100% galvanic isolation and Zoom +/-
USB with 100% galvanic isolation



Configurations



Configurations

